

## VISIBLE SPECTROPHOTOMETER

## Feature







## Parameters

| Model                          | 721  | 722   | 723N   | 723PC  |
|--------------------------------|--|---|--|--|
| Wavelength Range(nm)           | 340 ~ 1000   | 325 ~ 1000  | 325-1000   | 325-1000   |
| Transmittance Range(%)         | 0 ~ 199.9  | 0~199.9   | 0~125.0  | 0~125.0  |
| Absorbance Range A             | -0.301~1.999   | -0.301~1.999  | -0.097~2.5   | -0.097~2.5   |
| Scope fo Concentration Display | 0~1999   | 0~1999  | 0~1999   | 0~1999   |
| Wavelength Accuracy            | ±3nm   | ±2nm  | ±1nm   | ±1nm   |
| Spectrum Bandwidth             | 5nm  | 4nm   | 4nm  | 4nm  |
| Measure Precision              | 1%T  | 0.5%T   | 0,5%T  | 0,5%T  |
| Stray Light                    | ≤1%T (360nm)   | ≤0.3%T (220nm)  | ≤0.3%T (360nm)   | ≤0.3%T (360nm)   |
| Transmittance Accuracy         | ±1%  | ±0.5%   | ±0.5%  | ±0.5%  |
| Transmittance Repeatability    | 0.30%  | 0.20%   | 0.20%  | 0.20%  |
| Characteristics                | Unique and attactive appearance, has T.A.C.F function  | Unique and attactive appearance, Special LCD with backlight and display function.                           | Has the function of self regulating 0%T&100%T  | Has the function of self regulating 0%T&100%T  |
|                                | Bigger sample room suit for 10cm cell, specifically used for testing the low content solution by sugar plants, power plants, waterworks, chemical industry, etc. | With the asvantage of hight accuracy wavelength, nice monochromatic, lower stray light etc.                 | The patented unique scanning device can make high accurav of test reading and good repeatability and stability | The patented unique scanning device can make high accuracy of test reading and good repeatability and stability  |
|                                | Special LCD with backlight and display function.   | With PC measuring system, hight accuracy in T-A switching over, Concentration factor setting, concentration | Bigger sample room suit for 10cm cell  | Bigger sample room suit for 10cm cell  |
|                                | Has the function of self regulating 0%T&100%T  | Has the function of self regualating 0%T&100%T  | With the high grade function of concentration direct-reading, GOTO λ   | With the high grade function of concentration direct-reading, GOTO λ   |
|                                | "KEY F"-with the patented function of regulation   | "KEY F"-with the patented function of regualating   | Has printing function can be equipped with a printer to print at a fixed time                                  | Has printing function can be equipped with a printer to print at a fixed time  |
|                                |  | Can print data by serial data printer.  |  | With on-line software, can do quantitative analysis and qualitative analysis for sample through computer controlled  |
|                                |  | It can match 10cm cell holder<br>and 2cm,3cm,5cm and 10cm cell<br>to expend measuring range                 |  | With on-line software, can do quantitative analysis and qualitative analysis for sample through computer controlled instruments spectral scans can be done |

www.faithful.cc -