FAITHFUL® FAITHFUL INSTRUMENT (Hebei) CO.,Ltd

CM SERIES CLARIFY TESTER

• They test clarify of injection and bottled medicine liquid.

Applicable standards

- •National standard(Chinese codex)/2015
- •Health department standard

• Corporate standard(Clarify tester)

CM-1 CLARIFY TESTER

Features

- •The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- •Flash eliminated. Illumination is adjustable.
- •Time and alarm automatically. Test continuously.
- •Used by one person.

Specifications

- •Lamp three-primary-color fluorescent lamp
- •Illumination range $(1000\sim4000)$ LUX
- •Power 220V/ 50Hz/ 30W
- •Time range $(1\sim99)$ s



•Dimension (66*18*50) cm

CM-2 CLARIFY TESTER



Features

- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- •Flash eliminated. Illumination is adjustable.
- Time and alarm automatically. Test continuously.
- •Used by two person at the same time.

Specifications

- •Lamp three-primary-color fluorescent lamp
- •Illumination range (1000~4000) LUX
- •Time range $(1\sim99)$ s
- ●Power 220V/ 50Hz/ 50W
- •Dimension (66*38*50) cm