

## Pharmacy Refrigerator

### FSF-5V316

Pharmacy refrigerators are specially designed to store medicines, vaccines, reagents and biomedical products.



#### Features

- Adjustable temperature range 2-8°C
- Microprocessor control, digital display, and temperature adjustment with an increment of 0.1°C  
Large digital display screen
- Temperature variation within  $\pm 3^{\circ}\text{C}$
- Permanently lubricated cooling fan for safety and longevity
- High efficiency specialty compressor with known field reliability
- Forced- air cooling with optimized air distribution system designed

#### Safety

- Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 8 hours without AC power.
- Visual and audible alarm system: High/low temperature, thermostat failure alarm, power failure alarm, door ajar.
- Power failure protection: turn on delay of the cooling system when power failure, restart delay protection.
- Keyboard lock, password protection to avoid set parameter randomly.
- Wide available voltage range: 187V~242V.

#### Humanized Design

- Optimized space utilization design with multi - level and adjustable shelf height to accommodate storage requirement of pharmaceutical products of various kinds.
- Standard with Auto on/off LED Interior light with manual switch on control panel.
- Safety lock to prevent unauthorized access
- 4 casters with stopper and 2 adjustable feet for easy moving and fixation.

**Reliable control system Inside temperature 2~8°C**

**Forced-air cooling system Frost-free**

**Service-free drainage system**

**Wide voltage range 187V~242V available**

**5 visual and audible alarms**

<b>Performance</b>		<b>Refrigeration System</b>	
Temperature Range	+2~+8°C	Compressor	1
Ambient Temperature	10~32°C	Refrigerant	R290
Controller	Microprocessor	Cooling Type	Forced-air cooling
Display	Digital Display	Defrost Type	Frost free
<b>Electrical Data</b>		<b>Alarm</b>	
Power Supply(V/Hz)	220/50~60; 115/60	High/Low Temperature	Y
Rated Power(W)	240W	Power Failure	Y
Electrical Current(A)	1.2	Door ajar	Y
Power Consumption	1.78kWh/24h	Thermostat failure	Y
<b>Dimensions</b>		<b>Accessories</b>	
Capacity	316L	Caster	4
Interior Size(W*D*H)	540*425*1380m	Lamp	LED
Exterior Size(W*D*H)	660*555*1915m m	Door	1/Glass door
Package Size(W*D*H)	710*640*1970m	Shelf	5
N.t./G.t.(approx)	92.5/100kg	Test hole	1/25mm
Loading Qty(20'/40'/40H)	27/54/54	Temperature recorder	Optional