

Catalogue of Balance

FA-N Series Internal Calibration Analytical Balance



Characteristics:

The visual appearance of the product is blended into the simple concept, stylish and generous. Specially makes the fine pouring cast-aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capacity. The panel with five operational buttons, sliding glass windscreen provides big space, insightful vision and operates easily. We use high precision magnet sensor to keep the weighing more fast and stable.

Main Technical Parameter:

Model	FA1004N	FA1104N	FA1204N	FA2004N	FA2104N	FA2204N
Capacity(g)	100	110	120	200	210	220
Resolution(g)	0.0001					
Min. weighing(g)	0.0004					
Stable Time	≤3s					
Display	LCD(white back light with font)					
Pan Size	Ø80mm					
Operation Temp.	5-25℃					
Repeat Ability	±0.0002g					
Liner	±0.0003g					
Draft shield size	240mm*190mm*265mm					
Cal. weight	Internal Calibration					
Option	RS232C & RS485 、 Printer					
Packing Size	Inside:415*285*430mm			Outside:475*315*485mm		
G. weight	10kg					