

03 HMC-MINI-6KS/6K HMC-MINI Centrifuge

Features

1. Composite shell, shatter-resistant
 2. The electronic braking system for rapid sample processing
 3. Quick release rotor tub, allows for easy rotor exchange
 4. HMC-6KS Two speed options: 4,000 & 6,000rpm
 5. HMC-6KS: 3-in-1 rotor design, high capacity and time-saving
--6x1.5/2.0ml tubes, --6x0.5ml tubes, --2-8x0.2ml PCR strip tubes
- Optional rotor: 4x5ml Eppendorf tube



01 HMC-6KS

Patented 3-in-1 rotor design



02 HMC-6KS 3-in-1 rotor



HMC-6K

Specification

Model	HMC-6K	HMC-6KS
Rotor	6x1.5ml/2.0ml rotor 2x8x0.2ml strip rotor	3-in-1 rotor (*)
Speed	6,000rpm	4,000 & 6,000rpm
Max.RCF	2,000 x g	1,200 & 2,000 x g
Dimensions	175x148x118mm	135x127x120mm
Weight	0.7kg	0.7kg
Power	AC100-240V	AC100-240V

*3-in-rotor ---- 6x1.5/2.0ml, 6x0.5ml tubes and 2-8x0.2ml PCR strip rotor

Ordering information

No.	Code	Description
01	AS-08060-00	HMC-6K, AC100-240V, 50/60Hz
02	AS-08030-00	HMC-6KS, AC100-240V, 50/60Hz

Accessories

HMC-6K		HMC-6KS	
AS-08061-01	Allen wrench	AS-08031-01	6x1.5/2.0ml, 6x0.5ml tubes & 2-8x0.2ml PCR Strip rotor
AS-08061-02	6 x 1.5/2.0ml rotor		
AS-08061-03	2 x 8 x 0.2ml strip rotor		
AS-08061-04	Slide rotor		
AS-08061-05	Adapters for 0.5/0.6ml tubes		
AS-08061-06	Adapters for 0.2ml tubes		