

## FSF-SPCT Horizontal Electrophoresis

### Characteristic:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.



### Specifications:

Gel Size (W×L)	Sample volume	Buffer Volume	Dimension (L×W×H)	Net weight
60×60mm; 120×60mm; 60×120mm 120×120mm	3, 6,8,11,13, 18 ,25.	~550ml.	300×160×75 (mm)	1.1 (kg)

## FSF-SPDT Horizontal Electrophoresis Cell

### Characteristic:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.



### Specifications:

Gel Size (W×L)	Sample volume	Buffer Volume	Dimension (L×W×H)	Net weight
140×140mm; 140×70mm; 70×70mm	1, 2, 3, 6, 14, 8,17, 13, 29	~900ml	390×187×63 (mm)	2.0 (kg)